


Polypectomy Snares

Single Use



Manufactured by:

Medi-Globe

Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by:

Medi-Globe
The Spirit of Care.™

Polypectomy Snares

Single Use


Standard Snare

Features:

- 360° Rotating Versions
- Continuous Braided Loop Wire
- Single Use

Benefits:

- Easier Polyp Snaring
- Smooth Cut
- Safety and Convenience



Rotatable Snare

Order Information:

Rotatable Polypectomy Snares:

Catalog No.	Description	Snare Diameter	Tube Diameter	Length
Symmetrical				
GPS-01-15-180	Oval Loop	15 mm	2.5 mm	180 cm
GPS-01-15-230	Oval Loop	15 mm	2.5 mm	230 cm
GPS-11-15-180	Oval Loop, Insulated Tip	15 mm	2.5 mm	180 cm
GPS-11-15-230	Oval Loop, Insulated Tip	15 mm	2.5 mm	230 cm
GPS-01-25-180	Oval Loop	25 mm	2.5 mm	180 cm
GPS-01-25-230	Oval Loop	25 mm	2.5 mm	230 cm
GPS-11-25-180	Oval Loop, Insulated Tip	25 mm	2.5 mm	180 cm
GPS-11-25-230	Oval Loop, Insulated Tip	25 mm	2.5 mm	230 cm
GPS-01-35-180	Oval Loop	35 mm	2.5 mm	180 cm
GPS-01-35-230	Oval Loop	35 mm	2.5 mm	230 cm
GPS-11-35-180	Oval Loop, Insulated Tip	35 mm	2.5 mm	180 cm
GPS-11-35-230	Oval Loop, Insulated Tip	35 mm	2.5 mm	230 cm

Rotatable Polypectomy Snares:

Catalog No.	Description	Snare Diameter	Tube Diameter	Length
Hexagonal				
GPS-02-25-180	Hexagonal Loop	25 mm	2.5 mm	180 cm
GPS-02-25-230	Hexagonal Loop	25 mm	2.5 mm	230 cm
GPS-12-25-180	Hexagonal Loop, Insulated Tip	25 mm	2.5 mm	180 cm
GPS-12-25-230	Hexagonal Loop, Insulated Tip	25 mm	2.5 mm	230 cm

Polypectomy Snares:

Catalog No.	Description	Snare Diameter	Tube Diameter	Length
Symmetrical				
GPS-21-15-180	Oval Loop	15 mm	2.5 mm	180 cm
GPS-21-15-230	Oval Loop	15 mm	2.5 mm	230 cm
GPS-31-15-180	Oval Loop, Insulated Tip	15 mm	2.5 mm	180 cm
GPS-31-15-230	Oval Loop, Insulated Tip	15 mm	2.5 mm	230 cm
GPS-21-25-180	Oval Loop	25 mm	2.5 mm	180 cm
GPS-21-25-230	Oval Loop	25 mm	2.5 mm	230 cm
GPS-21-25-300	Oval Loop	25 mm	2.5 mm	300 cm
GPS-31-25-180	Oval Loop, Insulated Tip	25 mm	2.5 mm	180 cm
GPS-31-25-230	Oval Loop, Insulated Tip	25 mm	2.5 mm	230 cm
GPS-21-35-180	Oval Loop	35 mm	2.5 mm	180 cm
GPS-21-35-230	Oval Loop	35 mm	2.5 mm	230 cm
GPS-31-35-180	Oval Loop, Insulated Tip	35 mm	2.5 mm	180 cm
GPS-31-35-230	Oval Loop, Insulated Tip	35 mm	2.5 mm	230 cm
GPS-21-50-180	Oval Loop	50 mm	2.5 mm	180 cm
GPS-21-50-230	Oval Loop	50 mm	2.5 mm	230 cm

Polypectomy Snares:

Catalog No.	Description	Snare Diameter	Tube Diameter	Length
Hexagonal				
GPS-22-25-180	Hexagonal Loop	25 mm	2.5 mm	180 cm
GPS-22-25-230	Hexagonal Loop	25 mm	2.5 mm	230 cm
GPS-32-25-180	Hexagonal Loop, Insulated Tip	25 mm	2.5 mm	180 cm
GPS-32-25-230	Hexagonal Loop, Insulated Tip	25 mm	2.5 mm	230 cm
GPS-22-35-180	Hexagonal Loop	35 mm	2.5 mm	180 cm
GPS-22-35-230	Hexagonal Loop	35 mm	2.5 mm	230 cm

Medi-Globe handles are directly compatible with the Olympus and Microvasive electro-surgical generators. To ensure compatibility with other types of electro-surgical generators, please contact us. We recommend the use of compatible Medi-Globe active cords. You will find the corresponding Medi-Globe active cords under "Accessories".