

RotaSnare®

Polypectomy Snares - Reusable



Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by:

 **Medi-Globe**[®]
The Spirit of Care.[™]

RotaSnare®

Polypectomy Snares - Reusable

Features:

- 360° Rotating Loop
- Continuous Braided Loop Wire
- Complete Component Detachability
- Autoclavable

Benefits:

- Easier Polyp Snaring
- Smooth Cut and Long Lasting Loop
- Economical Component Replacement
- Easy and Effective Sterilisation

Order Information:

RotaSnare® Symmetrical Snares, rotatable


Catalog No.	Description	Tube Diameter	Snare Diameter	Length
Polypectomy Snares (without handle and tube):				
713360	Oval	2.5 mm	15 mm	180 cm
713460	Oval	2.5 mm	15 mm	230 cm
714360	Oval, with insulated tip	2.5 mm	15 mm	180 cm
714460	Oval, with insulated tip	2.5 mm	15 mm	230 cm
715360	Oval	2.5 mm	25 mm	180 cm
715460	Oval	2.5 mm	25 mm	230 cm
716360	Oval, with insulated tip	2.5 mm	25 mm	180 cm
716460	Oval, with insulated tip	2.5 mm	25 mm	230 cm
717360	Oval	2.5 mm	35 mm	180 cm
717460	Oval	2.5 mm	35 mm	230 cm
718360	Oval, with insulated tip	2.5 mm	35 mm	180 cm
718460	Oval, with insulated tip	2.5 mm	35 mm	230 cm
721360	Hexagonal	2.5 mm	25 mm	180 cm
721460	Hexagonal	2.5 mm	25 mm	230 cm
722360	Hexagonal, with insulated tip	2.5 mm	25 mm	180 cm
722460	Hexagonal, with insulated tip	2.5 mm	25 mm	230 cm

Asymmetrical Snares, non-rotatable


Catalog No.	Description	Tube Diameter	Snare Diameter	Length
Polypectomy Snares (without handle and tube):				
703360	Asymmetrical	2.5 mm	15 mm	180 cm
703460	Asymmetrical	2.5 mm	15 mm	230 cm
705360	Asymmetrical	1.8 mm	25 mm	180 cm
706360	Asymmetrical	2.5 mm	25 mm	180 cm
706460	Asymmetrical	2.5 mm	25 mm	230 cm
708360	Asymmetrical, with insulated tip	2.5 mm	25 mm	180 cm
708460	Asymmetrical, with insulated tip	2.5 mm	25 mm	230 cm
709360	Asymmetrical	2.5 mm	40 mm	180 cm
709460	Asymmetrical	2.5 mm	40 mm	230 cm
710360	Asymmetrical, with insulated tip	2.5 mm	40 mm	180 cm
710460	Asymmetrical, with insulated tip	2.5 mm	40 mm	230 cm



Oval



Hexagonal



Asymmetrical

Oval Snares, non-rotatable

Catalog No.	Description	Tube Diameter	Snare Diameter	Length
Polypectomy Snares (without handle and tube):				
723360	Onion shaped loop	2.5 mm	50 mm	180 cm
723460	Onion shaped loop	2.5 mm	50 mm	230 cm
726360	Oval, monofil	2.5 mm	25 mm	180 cm
726460	Oval, monofil	2.5 mm	25 mm	230 cm


Hexagonal Snares, non-rotatable

Catalog No.	Description	Tube Diameter	Snare Diameter	Length
Polypectomy Snares (without handle and tube):				
720360	Hexagonal	2.5 mm	25 mm	180 cm
720460	Hexagonal	2.5 mm	25 mm	230 cm
725360	Hexagonal	2.5 mm	35 mm	180 cm
725460	Hexagonal	2.5 mm	35 mm	230 cm

Related Accessories:

703363	Insertion Tube for Catalog no.: 703360, 706360, 707360, 708360, 709360, 710360, 711360, 712360
703463	Insertion Tube for Catalog no.: 703460, 706460, 707460, 708460, 709460, 710460, 711460, 712460
705363	Insertion Tube for Catalog no.: 705360
713363	Insertion Tube for Catalog no.: 713360, 714360, 715360, 716360, 717360, 718360, 720360, 721360, 722360, 723360, 725360 und 726360
713463	Insertion Tube for Catalog no.: 713460, 714460, 715460, 716460, 717460, 718460, 720460, 721460, 722460, 723460, 725460 und 726460
100700	Handle for all Polypectomy Snares (except for 709...and 710...)
100703	Handle for 709...and 710...

100700



100703



Snare Wire, Insertion Tube and Handle are required for a complete reusable unit.

Medi-Globe handles are compatible with Olympus and Microvasive electro-surgical generators. To ensure compatibility with other types of electro-surgical generators, please contact us or your local distributor.