

Aspiration Needles

Single Use



Aspiration Needles


for Bronchoscopy and Gastroscopy - Single Use

Features:

- Coaxial Luer-Lock for Vacuum Syringe
- Lock Mechanism
- Easy Grip Handle
- 15° Beveled Stainless Steel Needle
- Single Use

Benefits:

- Easy Aspiration of Cells
- Secure Needle Introduction and Projection
- Safe and Controllable Application
- Maximum Piercing Capability
- Safety and Convenience



Needle with Side Hole

Order Information:

Catalog No.	Description	Needle Diameter	Needle Length	Outer Diameter	Total Length
GAN-01-19-100	Aspiration Needle for Bronchoscopy	0.8 mm / 21 Gauge	12 mm	1.9 mm	100 cm
GAN-02-19-180	Aspiration Needle for Gastroscopy	0.8 mm / 21 Gauge	12 mm	1.9 mm	180 cm

Manufactured by:

Medi-Globe
 Medi-Globe GmbH
 Medi-Globe-Strasse 1-5
 83101 Achenmühle, Germany
 Phone: ++49 8032 973 379
 Fax: ++49 8032 973 339
 e-mail: sales@medi-globe.de
 www.medi-globe.de

Distributed by: