


Stone Extraction Baskets

Reusable



Stone Extraction Baskets

Reusable

Features:

- Continuous Braided Wire and Baskets
- Adjustable Sealing Cap
- Ergonomic Instrument Holder
- All Components are Autoclavable

Benefits:


- Minimizes Basket / Cable Breakage
- Contrast Compatibility
- Simple and Safe Manipulation of the Basket
- Easy and Effective Sterilization

Order Information:

Stone Extraction Basket Set – consists of Basket, Teflon® Tube, Instrument Holder and Sealing Cap

Catalog No.	Description	Outer Diameter	Basket Diameter
Basket with Tube 180 cm			
606420	4-Wire Basket	1.9 mm	20 mm
606425	4-Wire Basket	1.9 mm	25 mm
606430	4-Wire Basket	1.9 mm	30 mm
608620	6-Wire Basket	1.9 mm	20 mm
608625	6-Wire Basket	1.9 mm	25 mm
608630	6-Wire Basket	1.9 mm	30 mm
605420	4-Wire Basket	2.5 mm	20 mm
605425	4-Wire Basket	2.5 mm	25 mm
605430	4-Wire Basket	2.5 mm	30 mm
607620	6-Wire Basket	2.5 mm	20 mm
607625	6-Wire Basket	2.5 mm	25 mm
607630	6-Wire Basket	2.5 mm	30 mm
Basket with Tube 280 cm			
615420	4-Wire Basket	2.5 mm	20 mm
615425	4-Wire Basket	2.5 mm	25 mm
615430	4-Wire Basket	2.5 mm	30 mm
615620	6-Wire Basket	2.5 mm	20 mm
615625	6-Wire Basket	2.5 mm	25 mm

As an alternative to the instrument holder the handle order no. 100600 can be used.



Manufactured by:

Medi-Globe

Medi-Globe GmbH

Medi-Globe-Strasse 1-5

83101 Achenmühle, Germany

Phone: ++49 8032 973 379

Fax: ++49 8032 973 339

e-mail: sales@medi-globe.de

www.medi-globe.de

Distributed by: