




# E.R.C.P. Catheters

Single Use



# E.R.C.P. Catheters

Single Use


## Features:

- "Pre-Curved"
- Metal Tips/ Metal Rings
- Guide Wire Compatibility
- Color Markings
- Removable Stylet (1.2 m)
- Single Use

## Benefits:

- Easier Positioning to the Papilla
- Maximum X-Ray Visualization
- Contrast Injection while Guide Wire in Place
- Macroscopically Controlled Cannulation Depth
- Kink Free Cannulation
- Safety and Convenience

## Order Information:



### E.R.C.P. Catheter

Catalog No.	Description	Outer Diameter	Guide Wire
GER-01-18-035	Standard Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)
GER-02-18-035	5 mm Tapered Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)
GER-03-18-021	10 mm Tapered Tip "Pre-Curved"	1.8 mm	.021" (0.53 mm)
GER-04-18-035	Metal Cone Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)
GER-05-18-035	Round Metal Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)
GER-06-18-035	"Bullet Type" Metal Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)
GER-07-18-035	Metal Ball Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)

Manufactured by:



Medi-Globe GmbH  
 Medi-Globe-Str. 1-5  
 83101 Achenmühle, Germany  
 Phone: ++49 8032 973 379  
 Fax: ++49 8032 973 339  
 e-mail: sales@medi-globe.de  
 www.medi-globe.de

Distributed by: