





E.R.C.P. Catheters

Reusable



E.R.C.P. Catheters

Reusable

Features:

- Teflon® Tip with Memory-Effect
- Color Markings
- Metal Tips / Metal rings
- Guide Wire Compatibility
- Removable Stylet (1.2 m)
- Autoclavable

Benefits:

- Allows Pre-Shaping of the Tip
- Macroscopically Controlled Cannulation Depth
- Maximum X-Ray Visualization
- Contrast Injection while Guide Wire in Place
- Kink Free Cannulation
- Easy and Effective Sterilization

Order Information:



E.R.C.P. Catheter

Catalog No.	Description	Outer Diameter	Accepted Guide Wire
522420	Tapered Tip 10 mm	1.8 mm	.021" (0.53 mm)
522421	Cone Tapered Tip 5 mm	1.8 mm	.021" (0.53 mm)
524420	Cone Metal Tip	1.8 mm	.021" (0.53 mm)
526420	Round Metal Tip	1.8 mm	.021" (0.53 mm)
528420	Tapered Metal Tip	1.8 mm	.021" (0.53 mm)
520420	Standard	1.8 mm	.035" (0.89 mm)
521420	Standard	2.2 mm	.035" (0.89 mm)
523420	Tapered Tip 10 mm	2.2 mm	.035" (0.89 mm)
523421	Cone Tapered Tip 5 mm	2.2 mm	.035" (0.89 mm)
525420	Cone Metal Tip	2.2 mm	.035" (0.89 mm)
527420	Round Metal Tip	2.2 mm	.035" (0.89 mm)
529420	Tapered Metal Tip	2.2 mm	.035" (0.89 mm)
532420	Bullet-Tip	1.8 mm	.035" (0.89 mm)
533420	Metal Ball Tip	1.8 mm	.035" (0.89 mm)
534420	Minor E.R.C.P. Catheter, Metal Cannula Length 1 mm	1.8 mm	
535420	Tapered Tip Length 30 mm, Diameter 2.4 mm Tapered to 1.85 mm	2.4 / 1.85 mm	.035" (0.89 mm)

Manufactured by:

Medi-Globe

Medi-Globe GmbH
 Medi-Globe-Strasse 1-5
 83101 Achenmühle, Germany
 Phone: ++49 8032 973 379
 Fax: ++49 8032 973 339
 e-mail: sales@medi-globe.de
 www.medi-globe.de

Distributed by: