

Foreign Body Retrievers

Reusable



Manufactured by:

Medi-Globe

Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by:

Medi-Globe
The Spirit of Care.™

Foreign Body Retrievers

Reusable

Features:

- Wide Opening of Graspers -----> Large Variety of Applications
- Stationary Four Prong Graspers -----> Eliminates the Retraction of the Grasper when Closing (Catalog No. 802...&. 803....)
- Detachable Distal Tips (Forceps) or -----> Optimal Cleaning and Sterilization Detachable Handle and Grasper (Retrievers)
- Ground Distal Spiral -----> Easy Passage through Working Channel
- All Components are Autoclavable -----> Easy and Effective Sterilization (Except Catalog No. 802... & 803...: ETO gas only)

Benefits:

Catalog No.	Description	Outer Diameter	Length
820360	Alligator Forceps	1.8 mm	180 cm
820460	Alligator Forceps	1.8 mm	230 cm
821360	Alligator Forceps	2.3 mm	180 cm
821460	Alligator Forceps	2.3 mm	230 cm

Order Information:

Catalog No.	Description	Outer Diameter	Length
802360	Foreign Body Retrievers, 4-prong grasper with hooks, stationary	2.6 mm	180 cm
802460	Foreign Body Retrievers, 4-prong grasper with hooks, stationary	2.6 mm	230 cm
803360	Polyp Retrievers, 4-prong grasper with eyelets, stationary	2.6 mm	180 cm
803460	Polyp Retrievers, 4-prong grasper with eyelets, stationary	2.6 mm	230 cm
804360	Foreign Body Retrievers, 4-prong grasper with hooks	2.5 mm	180 cm
804460	Foreign Body Retrievers, 4-prong grasper with hooks	2.5 mm	230 cm
805360	Polyp Retrievers, 4-prong grasper with eyelets	2.5 mm	180 cm
805460	Polyp Retrievers, 4-prong grasper with eyelets	2.5 mm	230 cm
807360	Foreign Body Retrievers, forked jaws (2:2) grasper	2.5 mm	180 cm
807460	Foreign Body Retrievers, forked jaws (2:2) grasper	2.5 mm	230 cm
812360	Foreign Body Retrievers, forked jaws (1:2) grasper	2.5 mm	180 cm
812460	Foreign Body Retrievers, forked jaws (1:2) grasper	2.5 mm	230 cm
816360	Rat Tooth Forceps	1.8 mm	180 cm
816460	Rat Tooth Forceps	1.8 mm	230 cm
817360	Rat Tooth Forceps	2.3 mm	180 cm
817460	Rat Tooth Forceps	2.3 mm	230 cm

Shapes of Foreign Body Retrievers, Polyp Retrievers and Foreign Body Forceps

