

Guide Wires

Single Use



Guide Wires

Single Use

Features:

- Nitinol Versions
- Hydrophilic Coated Versions
- Teflon® Coated
- Two Color Markings and Radiopaque Tip
- Single Use or Autoclavable

Benefits:

- Kink Resistant
- Easy, Low Friction Manipulation
- Usable during Sphincterotomy
- Enhanced Visualization
- Safety and Convenience or Price Reduction by Multi-Use

Order Information:

Catalog No.	Description	Diameter	Length
Single Use	Teflon® coated, stainless steel		
GGW-01-21-400	Fixed Core	.021" (0.53 mm)	400 cm
GGW-01-21-450	Fixed Core	.021" (0.53 mm)	450 cm
GGW-01-35-400	Fixed Core	.035" (0.89 mm)	400 cm
GGW-01-35-450	Fixed Core	.035" (0.89 mm)	450 cm
GGW-01-47-400	Fixed Core	.047" (1.20 mm)	400 cm
GGW-01-47-450	Fixed Core	.047" (1.20 mm)	450 cm
GGW-01-52-400	Fixed Core	.052" (1.32 mm)	400 cm
GGW-01-52-450	Fixed Core	.052" (1.32 mm)	450 cm
GGW-02-35-400	Movable Core	.035" (0.89 mm)	400 cm
GGW-02-35-450	Movable Core	.035" (0.89 mm)	450 cm

Catalog No.	Description	Diameter	Length
Autoclavable	Hydrophilic Teflon® Coated Nitinol Wire, two-color markings		
GGW-03-21-400	"Soft" Tip, 50 mm tapered radiopaque Tip	.021" (0.53 mm)	400 cm
GGW-03-21-450	"Soft" Tip, 50 mm tapered radiopaque Tip	.021" (0.53 mm)	450 cm
GGW-03-35-400	"Soft" Tip, 50 mm tapered radiopaque Tip	.035" (0.89 mm)	400 cm
GGW-03-35-450	"Soft" Tip, 50 mm tapered radiopaque Tip	.035" (0.89 mm)	450 cm
GGW-04-35-400	"Standard" Tip, 50 mm tapered radiopaque Tip	.035" (0.89 mm)	400 cm
GGW-04-35-450	"Standard" Tip, 50 mm tapered radiopaque Tip	.035" (0.89 mm)	450 cm

Manufactured by:



Medi-Globe-Str. 1-5

83101 Achenmühle, Germany

Phone: ++49 8032 973 379

Fax: ++49 8032 973 339

e-mail: sales@medi-globe.de

www.medi-globe.de

Distributed by: