


Injectra

Injection Needles - Single Use



 **Medi-Globe**[®]
The Spirit of Care.[™]

Injectra

Injection Needles - Single Use

Features:

- Top Down Injection
- Lock Mechanism
- Easy Grip Handle
- 15° Beveled Stainless Steel Needle
- Single Use

Benefits:

- Easy Injection of Sclerosant
- Secure Needle Introduction and Projection
- Safe and Controllable Operation
- Maximum Piercing Capability
- Safety and Convenience

Order Information:

Injection Needle

Catalog No.	Needle Diameter	Needle Length	Catheter Diameter	Catheter Length
GSN-01-19-180	0.7 mm / 22 Gauge	4.0 mm	1.9 mm	180 cm
GSN-01-19-230	0.7 mm / 22 Gauge	4.0 mm	1.9 mm	230 cm
GSN-02-23-180	1.0 mm / 19 Gauge	4.0 mm	2.3 mm	180 cm
GSN-02-23-230	1.0 mm / 19 Gauge	4.0 mm	2.3 mm	230 cm
GSN-03-19-120	0.7 mm / 22 Gauge	6.0 mm	1.9 mm	120 cm
GSN-03-19-180	0.7 mm / 22 Gauge	6.0 mm	1.9 mm	180 cm
GSN-03-19-230	0.7 mm / 22 Gauge	6.0 mm	1.9 mm	230 cm
GSN-04-23-120	1.0 mm / 19 Gauge	6.0 mm	2.3 mm	120 cm
GSN-04-23-180	1.0 mm / 19 Gauge	6.0 mm	2.3 mm	180 cm
GSN-04-23-230	1.0 mm / 19 Gauge	6.0 mm	2.3 mm	230 cm
GSN-05-19-180	0.5 mm / 25 Gauge	4.0 mm	1.9 mm	180 cm
GSN-05-19-230	0.5 mm / 25 Gauge	4.0 mm	1.9 mm	230 cm
GSN-06-19-180	0.5 mm / 25 Gauge	6.0 mm	1.9 mm	180 cm
GSN-06-19-230	0.5 mm / 25 Gauge	6.0 mm	1.9 mm	230 cm

Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH
Medi-Globe-Str. 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by: