


Fibrin Needles

Single Use



Fibrin Needles

Single Use

Features:

- Double Lumen Coaxiale "Needle in Needle" System
- Removable Safety Clip
- 15° Beveled Needle
- Easy Grip Handle
- Single Use

Benefits:

- > Precise Confluence of Fibrin/Thrombin at Intended Position
- > Secure Needle within the Catheter
- > Effective Needle Puncture
- > Easy Handling
- > Safety and Convenience

Order Information:



Catalog No.	Needle Diameter	Needle Length	Catheter Diameter	Total Length
GFN-01-04-180	1.2 mm	4.0 mm	2.5 mm	180 cm
GFN-01-04-230	1.2 mm	4.0 mm	2.5 mm	230 cm
GFN-01-06-180	1.2 mm	6.0 mm	2.5 mm	180 cm
GFN-01-06-230	1.2 mm	6.0 mm	2.5 mm	230 cm
GFN-01-08-180	1.2 mm	8.0 mm	2.5 mm	180 cm
GFN-01-08-230	1.2 mm	8.0 mm	2.5 mm	230 cm

Manufactured by:

Medi-Globe

Medi-Globe GmbH

Medi-Globe-Strasse 1-5

83101 Achenmühle, Germany

Phone: ++49 8032 973 379

Fax: ++49 8032 973 339

e-mail: sales@medi-globe.de

www.medi-globe.de

Distributed by: