


Biopsy Forceps

Single Use



Biopsy Forceps

Single Use

Features:

- Short Rigid Forceps Tips
- Large Jaw Opening
- Single Use

Benefits:

- Easy Passage through Angled Working Channels
- Facilitates Tangential Biopsies
- Safety and Convenience

Order Information:

Catalog No.	Description	Outer Diameter	Length	Minimum Working Channel Diameter
GBF-01-22-180	Biopsy Forceps, oval cups	2.4 mm	180 cm	2.6 mm
GBF-01-22-230	Biopsy Forceps, oval cups	2.4 mm	230 cm	2.6 mm

Manufactured by:

 **Medi-Globe**

Medi-Globe GmbH
Medi-Globe-Strasse 1-5
83101 Achenmühle, Germany
Phone: ++49 8032 973 379
Fax: ++49 8032 973 339
e-mail: sales@medi-globe.de
www.medi-globe.de

Distributed by: